

Performance of Diagnostical SARS-CoV-2 Mutation Detection Kits (RUO)*

A Validation Study with a Unique Cohort

Murat İşbilen¹, Serpil Bayraktar¹, Beyza Rümeysa Geç¹, Elif Yıldırım¹, Murat Sayan²



AIM

Various mutations in different regions of SARS-CoV-2 genome have been reported since the beginning of the COVID-19 pandemic [1-3], thus resulting in occurrence of new SARS-CoV-2 variants with increased transmissibility, disease severity and immune escape. The Centers for Disease Control and Prevention (CDC) and the World Health Organization (WHO) classified variants as variant of concern (VOC) or variant of interest (VOI) based on several criteria, such as their potential to affect clinical outcome [4, 5].

Identification of the mutations that have significant effect on the behaviour of the virus is crucial for the global control of the Covid-19. Reverse transcription real-time PCR (RT-qPCR) can be used as an alternative to gold standard, which is next-generation sequencing (NGS). The aim of this study is to assess that Diagnostical SARS-CoV-2 Mutation Detection Kits are capable of detecting the mutations in circulating SARS-CoV-2 variants by performing *in silico* analysis with data from GISAID database [6] and validating *in silico* results *in vitro* by sequencing assessment.

METHODS

We performed *in silico* and *in vitro* analyses to reveal the validity of the assays in detecting mutations in SARS-CoV-2 variants. We downloaded sequences of 46,546 SARS-CoV-2 specimens for 24 different variants that were classified as VOC or VOI by CDC or WHO, from the GISAID database [6]. Primer and probe sequences of the assays were aligned to these sequences and checked if there is any mismatches between the oligo sequences and the sample sequences.

A sample genome was assumed to be detected by the assay, in case all primers and probe sequences align to sample genome without a mismatch. Performances of the assays were calculated as the percent samples detected by the primer/probe set of the assay.

Validation study was performed by a unique cohort of SARS-CoV-2 specimens collected in Kocaeli University Hospital, Turkey. RT-qPCR experiments were performed according to the instructions explained in the user manual of DIAGNOSTICAL SARS-CoV-2 Mutation Detection Kits. The samples were also sequenced by NGS for the whole S gene, including the receptor-binding domain.

RESULTS

1. *in silico* Results

In silico analyses of the DIAGNOSTICAL SARS-CoV-2 Mutation Detection Assays showed that the assays have around 99% *in silico* detection rate (Table 1).

Table 1: <i>in silico</i> detection performance of Diagnostical Mutation Detection Kits			
Mutation	# of Samples Detected	# of Samples not Detected	<i>in silico</i> Detection Rate
N501Y*	14828	22	99.9%
E484Q*	2410	5	99.8%
T478K*	12684	33	99.7%
P681R*	15016	68	99.5%
E484K*	13854	43	99.7%
K417N**	15848	161	99.0%
P681H*	5277	25	99.5%
L452R**	19137	58	99.8%
K417T*	5062	6	99.9%
HV69/70del*	2479	41	98.4%
T547K**	4266	3	99.9%

* Samples from Omicron variant were excluded, as there are Omicron-specific mutations in the target region

** Part of the DIAGNOSTICAL SARS-CoV-2 High-Cov Multiplex Mutation Detection Assay, which is under development.

in silico detection rate of the assays were calculated based on the detection of mutant probe for the positive identification of mutant specimens. Note that the working principal of the DIAGNOSTICAL SARS-CoV-2 Mutation Detection Assays is based on the dual probe approach (Figure 1) for each target mutation to detect both wild-type and mutant alleles.




Figure 1: Molecular design of the DIAGNOSTICAL SARS-CoV-2 Mutation Detection Assays.

Because the assays differentiate a single nucleotide difference within sample sequences, they are sensitive to mismatches between the oligos and the sample genome. Therefore, any other mutations in the regions targeted by the assays may affect the performance of the assays significantly. To eliminate this possibility, a holistic approach is used to design the oligos by evaluating the mutational profile of different SARS-CoV-2 variants at the same time. The oligos were aligned on to the reference SARS-CoV-2 genome together with the mutations present in any of the variants and designed such that there is no overlap between the mutations and the target regions (Figure 2).

We analyzed mutational profiles of the variants to see if the combined evaluation of the assays may distinguish different




Figure 2: Variant specific mutations in K417N target region. Mutations are represented by the blue boxes for the variants in the y-axis. Position of the mutations on the genome indicated by the lines at the bottom of the figure. Red, green and yellow colors on the genome represent the regions targeted by forward primer, probe and reverse primer, respectively.




Figure 3: Mutation profiles of VOC and VOI variants in GISAID database. Each square represents the presence of a mutation in the right y-axis in a variant shown in the x-axis. The prevalence of the mutations in the variants is shown by the shades of purple/red color. The more purple/red the square is, the more prevalent the mutation is. The genomic positions of the mutations are shown in the left y-axis. The vertical gray line represents the S gene in the reference genome. The reference genome used is MN908947.3 [7].

variants (Figure 3). Alpha, Beta and Gamma variants can be identified by N501Y, K417N, E484K and HV69/70Del assays, as Delta, Kappa and Epsilon variants can be identified by E484Q, E484K, L452R and K417N assays. The Omicron variant can be distinguished from the Delta variant by the means of K417N and L452R mutations. It is also possible to distinguish the Omicron subvariants, BA.1 and BA.2 by using the K417N and the HV69/70del mutations (Figure 3).

2. *in vitro* Validation Results

In line with the *in silico* results, DIAGNOSTICAL SARS-CoV-2 Mutation Detection Kits showed a great concordance with NGS in clinical samples. All the clinical samples were genotyped correctly for all the mutations, except for two samples that were undefined by N501Y and HV69/70del assays (Table 2). Those two samples did not display any amplification for any of the mutant and wild-type assays of N501Y and HV69/70del mutations, because there were other minor mutations in the target regions on the genome of those samples (data not shown).

Table 2: *in vitro* detection performance of Diagnostical Mutation Detection Kits in the unique cohort.

RT-qPCR:	Mutant in NGS			WT in NGS			Sens.	Spec.
	Mut	WT	Und*	Mut	WT	Und*		
N501Y	109	0	0	58	1	0	100%	98.31%
HV69/70 deletion	95	0	1	0	71	0	98.96%	100%
E484K	10	0	0	0	158	0	100%	100%
L452R	28	0	0	0	247	0	100%	100%
K417N	100	0	0	0	15	0	100%	100%
E484Q	2	0	0	0	166	0	100%	100%
T547K	82	0	0	0	16	0	100%	100%

Mut: Mutant, WT: Wild-type, Und: Undefined, Sens: Sensitivity, Spec: Specificity

*Inconclusive due to another mutation in the target region

We performed another validation study for the new DIAGNOSTICAL High-Cov Multiplex Mutation Detection Assay that targets T547K, K417N and L452R mutations together with a SARS-CoV-2-specific Orflab sequence and a human endogenous control in five channels. The study showed that the multiplex assay could detect the mutations in the Omicron clinical samples (Table 3A) and the presence of the virus (Table 3B). The assay can be used to differentiate between BA.1 and BA.2 Omicron subvariants, as T547K mutation is BA.1 specific (Table 3C). Another validation study with the Delta samples is to be performed to assess the capability of the assay to distinguish the Delta and the Omicron variants.

Table 3: *in vitro* detection performance of Diagnostical SARS-CoV-2 High-Cov Multiplex Mutation Detection Assay in the unique cohort.

RT-qPCR:		Mutant in NGS		WT in NGS		Sens.	Spec.
		Mut	NoA	Mut	NoA		
T547K	56	0	0	10	0	100%	100%
K417N	65	0	0	0	0	100%	NA
L452R	0	0	0	65	0	NA	100%

RT-qPCR Sample (Orflab) s		NextStrain Clade (based on NGS sequences)			
		21K (BA.1)		21L (BA.2)	
Det	66				
NDet	0				
%Det	100%				
		55	0	1	10

Mut: Mutant, WT: Wild-type, NoA: Not Amplified, Und: Undefined, Sens: Sensitivity, Spec: Specificity, Det: Detected, NDet: Not Detected, %Det: Percent Detected.

CONCLUSIONS

The performance of DIAGNOSTICAL SARS-CoV-2 Mutation Detection Assays were assessed with the data collected from GISAID database and validated with a clinical sample set. Although the kits have high performance in detecting targeted mutations, other mutations in the target region may prevent effective detection of mutations. The new DIAGNOSTICAL SARS-CoV-2 High-Cov Multiplex Mutation Detection Assay is capable of distinguishing subvariants of the Omicron. These results suggest that real-time qPCR is an efficient method for the detection of SARS-CoV-2 mutations despite its limitations.

REFERENCES

- [1] Rambaut, A. et al. "Preliminary Genomic Characterisation of an Emergent SARS-CoV-2 Lineage in the UK Defined by a Novel Set of Spike Mutations." *Virological*, 18 Dec. 2020, virological.org/t/preliminary-genomic-characterisation-of-an-emergent-sars-cov-2-lineage-in-the-uk-defined-by-a-novel-set-of-spike-mutations/563.
- [2] Tegally, H. et al. Emergence and rapid spread of a new severe acute respiratory syndrome-related coronavirus 2 (sars-cov-2) lineage with multiple spike mutations in south africa. *medRxiv* (2020).
- [3] Voloch, C. M. et al. Genomic characterization of a novel sars-cov-2 lineage from rio de janeiro, brazil. *medRxiv* (2020).
- [4] Centers for Disease Control and Prevention. SARS-CoV-2 Variant Classifications and Definitions. <https://www.cdc.gov/coronavirus/2019-ncov/variants/variant-info.html>, Accessed March 28, 2022.
- [5] World Health Organization. Tracking SARS-CoV-2 variants. Retrieved from <https://www.who.int/en/activities/tracking-SARS-CoV-2-variants/>, Accessed March 28, 2022.
- [6] Bogner P, Capua I, Cox NJ, Lipman DJ, others A global initiative on sharing avian flu data. *Nature*. 2006;442(7106):981.
- [7] Nucleotide [Internet]. Bethesda (MD): National Library of Medicine (US); National Center for Biotechnology Information; [1988] -. Accession No. MN908947.3 Severe acute respiratory syndrome coronavirus 2 isolate Wuhan-Hu-1, complete genome, [cited 2021 Jan 05]. Available from: <https://www.ncbi.nlm.nih.gov/nucleotide/MN908947.3?report=fasta>