Instructions for Use

REF 05032025

Intended Use

Cirrus COVID-19 IgM/IgG Antibody Test is a rapid and qualitative immunochromatographic in vitro assay for differential detection of IgM and IgG antibodies to SARS-CoV-2 in human serum, plasma and whole blood samples.

The device is designed to aid in the determination of recent or previous exposure to SARS-CoV-2 virus, tracking the body's immunity status after infection by the SARS-CoV-2 virus. It only provides a preliminary result. Negative results do not exclude SARS-CoV-2 infection. Patient management decisions should combine negative results with clinical observations, patient history, and epidemiological information.

A positive result does not necessarily mean a current infection but represents a different stage of the disease after infection. Current infection should be confirmed by Real-Time Reverse Transcriptase (RT-PCR) or viral gene sequencing.

The test is *in vitro* diagnostic (IVD) device for professional use only.

Background

Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) emerged in early 2020 with human cases in Wuhan, China (Zhou et al., 2020; Zhu et al., 2020). It has rapidly rampaged worldwide, causing a pandemic of coronavirus disease (COVID-19) that ranges from fever and breathing difficulty to acute respiratory distress and death (Huang et al., 2020; Zhu et al., 2020) due to massive alveolar damage and progressive respiratory failure. On March 11, 2020, the World Health Organization (WHO) has declared the global outbreak of COVID-19 a pandemic associated with substantial morbidity and mortality. SARS-CoV-2 causes the most severe disease in older patients and people with co-morbidities, including heart disease, diabetes, and other health conditions (Wu and McGoogan, 2020).

Before 2019, six alpha- and beta-coronaviruses were known to cause respiratory diseases of different severity, including four common cold coronaviruses (229E, NL63, OC43, and HKU1) and two highly pathogenic coronaviruses (severe acute respiratory syndrome [SARS-CoV] and Middle East respiratory syndrome [MERS-CoV]), which emerged in 2003 and since 2012, respectively (As-siri et al., 2013; Huang et al., 2020). Importantly, with massive hospitalization rates and high mortality, SARS-CoV-2 remains a major threat to humankind and intervention strategies are being rapidly pursued.

Principle of the Test

Cirrus COVID-19 IgM/IgG Antibody Test is an antibody-capture immunochromatographic assay for the simultaneous detection and differentiation of IgM and IgG antibodies to SARS-CoV-2 virus. When the sample (human serum, plasma or whole blood) is added in the sample well, the specific antibodies, if any in the sample, will react with the SARS-CoV-2 specific antigens conjugated to colloidal gold. As the immunocomplex migrates along the strip, IgM and IgG antibodies are captured on the respective areas where anti-human IgM and IgG antibodies were immobilized and red/violet lines appear, indicating positive results. A colored procedure control line always forms at the end of the test area marked "C". If SARS-CoV-2 antibodies are absent in the sample, no red/violet line will form in the test lines (G and M), indicating a negative result.

Materials Provided

- 25 test devices, individually packed
- 1 dropper bottle containing buffer reagent (3 ml)
- 25 sample transfer pipette, in test packages
- Instructions for use

Other Materials Required

- Alcohol swab
- Lancets
- Timer

Storage and Stability

Store the kit components between 2°C-30°C (36°F-86°F). Do not freeze.

Keep away from direct sunlight, moisture and heat.

The kit can be used until the expiry date indicated on the label.

Warnings and Precautions

- Read the instructions carefully before using the test.
- This test is designed for in vitro diagnostic use and professional use only.
- Do not use beyond the expiry date indicated on the package.
- The test is for single use only. Do not reuse the test.
- Do not use the test if the protective package is damaged.
- Samples may be infectious. Do not pipette by mouth, wear disposable gloves and lab coat. Avoid splashing or aerosol formation. Wipe spills thoroughly using an appropriate disinfectant.
- Dispose all samples, kit components and other test materials in a biohazard container according to local regulations.
- A definitive clinical diagnosis should not be based on the results of a single test but should only be made by the physician after all clinical and laboratory findings have been evaluated.
- If any result is found doubtful then one should seek for confirmation using other test methods.
- The test should be performed between 15°C-30°C (59°F-86°F). If stored refrigerated, ensure that kit components are brought to operating temperature before testing.

Specimen Collection and Handling

- Human serum, plasma and whole blood specimens can be used as test samples.
- Collect the specimens aseptically by venipuncture or fingerstick ways.
- Follow your facility's routine procedures to collect samples.
- Whole blood samples should be tested immediately. Serum or plasma samples can be stored at 2°C-8°C upto 7 days or at -20°C upto 6 months before being tested. Avoid freeze and thaw cycles.
- If serum/plasma samples have to be transported for testing they should be in leak-free tubes and should be transported in cold chain (2°C-8°C).
- Only clear and clean specimens should be tested. Hemolyzed, contaminated, grossly lipemic or turbid samples should not be used as test material.
- Handle all specimens as they are infectious. Follow regulations and take necessary precautions when handling infectious materials.
- Discard all specimens and materials which come in contact with the specimens as infectious waste.

Test Procedure

1. Remove the testing device from the sealed pouch and place it on a flat, dry surface. The opened test kit should be better used within 4 hours.
2. For fingerstick or venous whole blood: Using a micropipette, collect the fingerstick whole blood or whole blood (20 µl). For serum/plasma: Using a pipette, collect the serum/plasma (10 µl). See the dots on the tip of the plastic pipette supplied with the test; the volume drawn upto the first dot is 10 µl, and upto the second dot is 20 µl.
3. Add the collected sample to upper area (close to test window) of sample well (marked S) on the test device without air bubbles.
4. Wait for a few seconds; add 2 drops (around 90 µl) of the sample buffer to the sample well (marked S) of the test device.
5. Read the results in 15-30 minutes. Strong positive specimens may produce positive result in the first minutes. Do not interpret results older than 30 minutes.

Quality Control

The device contains an internal quality control, a red/violet colored band in the control region (C). An external control is not provided with the test kit, however, good laboratory practice recommends the daily use of an external control to ensure proper device performance. External quality control samples should be tested according to the standard quality control requirements established by individual laboratory.

Interpretation of Results

Read the results in 15-30 minutes. Do not interpret the results older than 30 minutes!

Negative Result

A colored line only on the control zone (C) of the test area.

A negative result means there is no IgG and IgM antibodies against SARS-CoV-2 in the specimen or the antibody concentration is below the detection limit of the test.

Positive Results

A colored line on the G zone and a colored line on the C zone of the test area.

This result indicates that the tested sample contains IgG antibodies against SARS-CoV-2.

A colored line on the M zone and a colored line on the C zone of the test area.


This result indicates that the tested sample contains IgM antibodies against SARS-CoV-2.

A colored line on the G zone, a colored line on the M zone and a colored line on the C zone of the test area.

This result indicates that the tested sample contains IgG and IgM antibodies against SARS-CoV-2.

Invalid Results

When there is no colored line on the control "C" zone of the test area other lines are disregarded and the result is invalid. The sample should be retested on a new device or another sample should be collected.

Limitations of the Test

- The test qualitatively detects the presence of human IgG and IgM antibodies against SARS-CoV-2 in the specimen.
- A negative test result does not rule out the possibility of a SARS-CoV-2 infection!
- The reagent can only be used to determine the immune status of the body to SARS-CoV-2 after infection, but not directly to diagnose current SARS-CoV-2 infection!
- Although the test demonstrates superior accuracy in detecting antibodies against SARS-CoV-2, a low incidence of false results can occur. Therefore, other clinically available tests are required in case of questionable results. A definitive clinical diagnosis should not be based on the results of a single test but should only be made by the physician after all clinical and laboratory findings have been evaluated.

Performance of the Test

Sensitivity and Specificity

87 K2EDTA plasma samples of PCR confirmed Covid-19 patients were tested using CE marked and FDA EUA granted ELISA kits and their antibody status were determined. The samples then tested with Cirrus COVID-19 IgM/IgG Antibody Test and the sensitivity of the test was calculated as follows:

IgM (NCP) Sensitivity: 92.85 %

IgG (NCP) Sensitivity: 92.85 %

IgG (spike) Sensitivity: 95.23 %

The specificity of Cirrus IgM/IgG Antibody Test was evaluated with 126 plasma samples which had been collected before December 2019. Only one sample produced an equivocal IgG signal. Thus, the specificity of Cirrus IgM/IgG Antibody Test was calculated as:

Specificity: 99.2%

Repeatability and Reproducibility

Tests showed positive results with all positive samples and showed negative results with negative samples. There was no significant difference observed for the same sample when repeatedly testing 10 tests in the same batch.

No appreciable intra and inter lot variation were observed among different tests for each lot, different lots, different operators at different test sites in different time for the same sample.

The results demonstrated that the repeatability and reproducibility of Cirrus COVID-19 IgM/IgG Antibody Test are satisfactory.

References

1. Li, etc., Early Transmission Dynamics in Wuhan, China of Novel Coronavirus-Infected Pneumonia, DOI:10.1056/NEJMoa2001316.
2. Zhou, P., Yang, X.L., Wang, X.G., Hu, B., Zhang, L., Zhang, W., Si, H.R., Zhu, Y., Li, B., Huang, C.L., et al. (2020). A pneumonia outbreak associated with a new coronavirus of probable bat origin. Nature 579, 270–273. Zhu, N., Zhang, D., Wang, W., Li, X., Yang, B., Song, J., Zhao, X., Huang, B.,
3. Shi, W., Lu, R., et al.; China Novel Coronavirus Investigating and Research Team (2020). A Novel Coronavirus from Patients with Pneumonia in China, 2019. N. Engl. J. Med. 382, 727–733.
4. Huang, C., Wang, Y., Li, X., Ren, L., Zhao, J., Hu, Y., Zhang, L., Fan, G., Xu, J., Gu, X., et al. (2020). Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet 395, 497–506.
5. Wu, Z., and McGoogan, J.M. (2020). Characteristics of and Important Lessons From the Coronavirus Disease 2019 (COVID-19) Outbreak in China: Summary of a Report of 72314 Cases From the Chinese Center for Disease Control and Prevention. JAMA. https://doi.org/10.1001/jama.2020.2648.
6. Assiri, A., Al-Tawfiq, J.A., Al-Rabeeah, A.A., Al-Rabiah, F.A., Al-Hajjar, S., Al-Barrak, A., Flemban, H., Al-Nassir, W.N., Balkhy, H.H., Al-Hakeem, R.F. et al. (2013). Epidemiological, demographic, and clinical characteristics of 47 cases of Middle East respiratory syndrome coronavirus disease from Saudi Arabia: a descriptive study. Lancet Infect. Dis. 13, 752–761.

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Referanslar

- Li, et al., Early Transmission Dynamics in Wuhan, China of Novel Coronavirus-Infected Pneumonia, DOI:10.1056/NEJMoa2001316.
- Zhou, P., Yang, X.L., Wang, X.G., Hu, B., Zhang, L., Zhang, W., Si, H.R., Zhu, Y., Li, B., Huang, C.L., et al. (2020). A pneumonia outbreak associated with a new coronavirus of probable bat origin. Nature 579, 270–273. Zhu, N., Zhang, D., Wang, W., Li, X., Yang, B., Song, J., Zhao, X., Huang, B.,
- Shi, W., Lu, R., et al.; China Novel Coronavirus Investigating and Research Team (2020). A Novel Coronavirus from Patients with Pneumonia in China, 2019. N. Engl. J. Med. 382, 727–733.
- Huang, C., Wang, Y., Li, X., Ren, L., Zhao, J., Hu, Y., Zhang, L., Fan, G., Xu, J., Gu, X., et al. (2020). Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet 395, 497–506.
- Wu, Z., and McGoogan, J.M. (2020). Characteristics of and Important Lessons From the Coronavirus Disease 2019 (COVID-19) Outbreak in China: Summary of a Report of 72314 Cases From the Chinese Center for Disease Control and Prevention. JAMA. https://doi.org/10.1001/jama.2020.2648.
- Assiri, A., Al-Tawfiq, J.A., Al-Rabeah, A.A., Al-Rabiah, F.A., Al-Hajjar, S., Al-Barrak, A., Flimban, H., Al-Nassir, W.N., Balkhy, H.H., Al-Hakeem, R.F. et al. (2013). Epidemiological, demographic, and clinical characteristics of 47 cases of Middle East respiratory syndrome coronavirus disease from Saudi Arabia: a descriptive study. Lancet Infect. Dis. 13, 752–761.

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