### INSTRUCTION FOR USE

#### RapidFor™ SARS-CoV-2 Rapid IgG/IgM Test Kit Catalog Number: VSCD01-25

For detection of IgG and IgM antibodies against SARS-CoV-2  
Detection in Whole Blood/Serum/Plasma

### FOR IN VITRO DIAGNOSTIC USE

This instruction for use (IFU) must be read carefully prior to use. Instruction for use must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions for use.

### INTENDED USE

**RapidFor™ SARS-CoV-2 Itop/d IgG/IgM test** is a colloidal gold enhanced, rapid immunoassay for the qualitative detection of SARS-CoV-2 IgG and IgM antibody in serum, plasma and whole blood samples from patients with suspected SARS-CoV-2 infection. This product is for in vitro emergency use only during the pneumonia epidemic of the SARS-CoV-2 infection since December 2019. The test kit should be complied with the relevant requirements of the "Pneumonitis Diagnosis and Treatment Scheme for SARS-CoV-2" and "Pneumonitis Prevention and Control Scheme for SARS-CoV-2" in use. The SARS-CoV-2 is a new type of coronavirus belonging to the genus  $\beta$ . It has a capsule, and its particles are round or oval, often polymorphous, with a diameter of 60-140 nm. Its genetic characteristics are significantly different from SARS-CoV and MERS-CoV. Current research shows that it has more than 85% homology with bat SARS-like coronavirus (bat-SL-CoVZC45). In vitro isolation and culture, SARS-CoV-2 can be found in human respiratory epithelial cells in about 96 hours, while it takes about 6 days to isolate and culture in Vero E6 and Huh-7 cell lines.

### PRINCIPLE OF THE PROCEDURE

The SARS-CoV-2 Rapid IgG/IgM Test is a Colloidal gold enhanced capture immunoassay for the qualitative determination of SARS-CoV-2 IgG and IgM antibody in human serum, plasma, and whole blood samples. The product is pre-embedded with recombinant SARS-CoV-2 antigen and rabbit antibody on the nitrocellulose membrane, and coated with anti-human IgM, anti-human IgG and goat anti-rabbit antibody respectively on the T1 line, T2 line and C line on the nitrocellulose membrane. When a positive sample is tested, the SARS-CoV-2 IgM and/or IgG antibody combines with the colloidal gold labeled SARS-CoV-2 antigen to form a complex. Under the action of chromatography, the complex flows along the membrane, when it passes T1 line and/or T2 line it combines with the anti-human IgM and/or IgG to form a colloidal gold complex to show color, and the gold colloidal labeled rabbit antibody combines with Goat anti-rabbit antibody at C line to show color. Negative samples only show color on C line.

### REAGENTS AND MATERIALS SUPPLIED

- Test devices with desiccant in individual foil pouches
- 1 Sample diluent (4 ml/vial)
- 25 Disposable Droppers
- 25 Lancet
- 1 Instructions for use

### STORAGE AND STABILITY

- Storage: Store in a dry place at 2-30°C
- Shelf life: 24 months.
- The test device should be used as soon as possible after being removed from the aluminum foil bag.
- MFD date and EXP dates are marked on aluminum foil bag.


### SPECIMEN REQUIREMENTS

1. It is suitable for serum, plasma, or whole blood samples. The commonly used anticoagulants (heparin, EDTA or sodium citrate) have no effect on the results of this kit.
2. Samples should be collected according to routine clinical methods and avoid hemolysis.
3. If serum or plasma specimens are to be tested in 5 days, they should be refrigerated at 2°C-8°C. Storage at -20°C should not exceed 3 months. For long time storage, they should be 70°C cryopreservation and avoid repeated freezing and thawing (no more than 3 times).
4. Whole blood samples can be stored at 2°C-8°C kept in refrigerator for 3 days with no cryopreservation.
5. Restore the sample to room temperature before test.
6. Obvious hemolysis, lipohemia and jaundice samples should not be used.

7. Transportation of the samples should be sealed ice cup with ice or sealed foam box with ice.

### TEST PROCEDURE

1. Open the aluminum foil bag, take out the test card, and lay it flat on a flat operating table and number it.
2. Add 1 drop (10  $\mu$ l) of specimen (whole blood/serum/plasma) to the sample well of the test card using the plastic dropper provided according to the figure.
3. Then add 2 drops (80  $\mu$ l) of sample diluent to the sample well.

### LIMITATIONS

1. The test results of this product cannot be used as a basis for diagnosis. Comprehensive judgment should be made in combination with clinical symptoms, epidemiological conditions, and further clinical data.
2. In the early stage of infection, the test result may be negative because the SARS-CoV-2 antibody or low antibody level has not yet appeared in the sample.
3. Due to the limitation of the detection method, the negative test result of this reagent cannot exclude the possibility of infection.
4. This reagent can only qualitatively detect SARS-CoV-2 antibodies in human serum, plasma, and whole blood samples. It cannot determine the certain antibody content in the samples.

### PERFORMANCE CHARACTERISTIC

**RapidFor™ SARS-CoV-2 Rapid IgG/IgM Test Kit** performance characteristics have been determined by using clinical samples. The PCR results have been used for comparison and the data which is shown below is determined.

SARS-CoV-2 IgM	ELISA RESULTS		Total
	+	-	
RapidFor™ SARS-CoV-2 Rapid IgG/IgM Test	277	23	300
Total	282	339	611

Sensitivity;  $277 / (277+23) = \%92.3$ ,  
Specificity;  $306 / (306+5) = \%98.30$ ,  
Total coincidence rate;  $(277+306) / 611 = \%95.40$ ,

SARS-CoV-2 IgG	ELISA RESULTS		Total
	+	-	
RapidFor™ SARS-CoV-2 Rapid IgG/IgM Test	290	10	300
Total	296	317	613

Sensitivity;  $290 / (290+10) = \%96.6$ ,  
Specificity;  $307 / (307+6) = \%98.08$ ,  
Total coincidence rate;  $(290+307) / 613 = \%97.30$

### Cross-Reactivity

Cross-reactivity of the RapidFor™ SARS-CoV-2 IgG/IgM Test was evaluated using serum samples which contain antibodies to the pathogens listed below. 24 different pathogen was tested and no false positives were found with the following

Sample Categories	Sample Number
Influenza A virus IgG	3
Influenza B virus IgG	3
Respiratory syncytial virus IgG	3
Adenovirus IgG	3
Rhinovirus IgG	3
Human metapneumovirus IgG	3
Mycoplasma pneumoniae IgG	3
Chlamydia pneumoniae IgG	3
HCV IgG	3
Haemophilus influenzae IgG	3
HBV core antibody IgG	3
Bacterial pneumonia	3

Influenza A virus IgM	3
Influenza B virus IgM	3
Respiratory syncytial virus IgM	3
Adenovirus IgM	3
Rhinovirus IgM	3
Human metapneumovirus IgM	3
Mycoplasma pneumoniae IgM	3
Chlamydia pneumoniae IgM	3
HCV IgM	3
Haemophilus influenzae IgM	3
HBV core antibody IgM	3
Antinuclear antibodies (ANA)	3

### Interference;

SARS-CoV-2 antibody positive serum samples and SARS-CoV-2 antibody negative serum samples investigated with one of the following substances to added with specified concentrations in multiple replicates. No false Positives or false negatives were found with the following:

Name of Substance	Concentration
Ascorbic acid	40 mg/dL
Hemoglobin	1200 mg/dL
Bilirubin	10 mg/dL
Albumin	2000 mg/dL
Triglyceride	1000 mg/dL

### REFERENCES

- [1] Bedford J, Enria D, Giesecke J, et al. COVID-19: towards controlling of a pandemic[J]. Lancet,2020.
- [2] Andersen K G, Rambaut A, Lipkin W I, et al. The proximal origin of SARS-CoV-2[J]. Nature Medicine,2020.
- [3] Guan W, Ni Z, Hu Y, et al. Clinical Characteristics of Coronavirus Disease 2019 in China[J]. The New England journal of medicine,2020.
- [4] Li Z, Yi Y, Luo X, et al. Development and Clinical Application of a Rapid IgM IgG Combined Antibody Test for SARS-CoV-2 Infection Diagnosis[J]. Journal of Medical Virology,2020.
- [5] Xie X, Zhong Z, Zhao W, et al. Chest CT for Typical 2019-nCoV Pneumonia: Relationship to Negative RT-PCR Testing[J]. Radiology,2020:200343.

### INDEX OF SYMBOL

	Manufacturer
	Consult instruction for use
	Attention, see instruction for use
	In vitro diagnostic medical device
	Number of test
	Catalog number
	For single use only
	Storage temperature
	Lot number
	Expiry date

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## REFERANSLAR

- [1]Bedford J, Enria D, Giesecke J, et al. COVID-19: towards controlling of a pandemic[J]. Lancet,2020.
- [2]Andersen K G, Rambaut A, Lipkin W I, et al. The proximal origin of SARS-CoV-2[J]. Nature Medicine,2020.
- [3]Guan W, Ni Z, Hu Y, et al. Clinical Characteristics of Coronavirus Disease 2019 in China[J]. The New England journal of medicine,2020. [4] Li Z, Yi Y, Luo X, et al. Development and Clinical Application of a Rapid IgM IgG Combined Antibody Test for SARS-CoV-2 Infection Diagnosis[J]. Journal of Medical Virology,2020.
- [5]Xie X, Zhong Z, Zhao W, et al. Chest CT for Typical 2019-nCoV Pneumonia: Relationship to Negative RT-PCR Testing[J]. Radiology,2020:200343.

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